



ADVANCED MOBILITY SYSTEMS

621 Justus Drive · Kingston, Ontario · K7M 4H5
800.661.6716 tel·613.384.7460 fax·613.634.3990
www.amstilt.com email·info@amstilt.com



PowerTilt2™

PowerTilt2™

The PowerTilt2 is designed to provide independence in tilt positioning. It offers a balance between comfort, maneuverability, speed and power.

FEATURES

- Maximum user weight of 136kg – 300lbs
- Tilt range of 45 degrees starting from 3 degrees anterior tilt
- Full tilt reached within 10 seconds
- Minimum number of full tilts and returns on one charge – 100 with a 300lb user
- Incremental weight of less than 6Kg – 13.5 lbs
- May be operated by either the care giver or the care recipient
- Can be fully recharged within six hours

SAFETY

- Does not negatively affect stability
- Battery pack fused for critical system conditions
- Automatic charger to prevent overcharging or overheating
- Audible warning when battery is low
- Deep discharge prevention
- Full functionality during charging

BENEFITS

- Client has independence to change their tilt-in-space position
- Hand control can be positioned anywhere for ease of access
- Minimal weight increase as compared to the iTilt2
- Can be operated by caregiver, relieving stress and strains associated with tilt-in-space wheelchairs
- The iTilt2 maintains it's full range of adjustability
- Low-profile design

FEATURES

- May be ordered on all new iTilt2 Chairs
- Older iTilt2 chairs can be retrofitted at the factory

MADE IN CANADA

ONE CHAIR, INFINITE POSSIBILITIES ▶▶▶